



Bionano Solve v3.7.1 Release Notes

Document Number: 30456

Document Revision: B

Table of Contents

Legal Notice.....	3
Revision History.....	4
Introduction.....	4
Compatibility	4
Use of IrysSolve servers.....	4
References	5
Bionano EnFocus™ Fragile X Analysis Pipeline.....	5
Bionano EnFocus™ FSHD Analysis Pipeline	5
AOH/LOH Detection	5
Variant Allele Fraction estimation	5
Copy number analysis pipeline.....	5
<i>De novo</i> assembly pipeline.....	6
Hybrid scaffolding	6
Molecule Quality Report (MQR).....	6
Rare Variant Pipeline (RVP).....	6
SV confidence	6
Variant Annotation Pipeline (VAP).....	6
VCF conversion	7
Miscellaneous.....	7
Other known issues and limitations	7
Technical Assistance.....	9

Legal Notice

For Research Use Only. Not for use in diagnostic procedures.

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics Customer Support for the latest information.

BIONANO GENOMICS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Patents

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

Trademarks

The Bionano Genomics logo and names of Bionano Genomics products or services are registered trademarks or trademarks owned by Bionano Genomics in the United States and certain other countries.

Bionano Genomics®, Irys®, IrysView®, IrysChip®, IrysPrep®, IrysSolve®, Saphyr®, Saphyr Chip®, Bionano Access®, and Bionano EnFocus™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano Genomics is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano Genomics. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2022 Bionano Genomics, Inc. All rights reserved.

Revision History

Revision	Effective Date	Notes
A	11/03/2021	Initial release of document.
B	05/09/2022	Updated with additions for Solve 3.7.1 release

Introduction

This document describes the release of Bionano Solve 3.7. We provide an overview of the fixes and improvements of the Bionano Solve analysis tools and pipelines so that you may better understand the impact of moving to this version of our software. Should you have any questions, please contact support@bionanogenomics.com.

Bionano Tools and Bionano Solve are now combined and branded as Bionano Solve. Bionano Solve is installed on Saphyr Compute, Bionano Compute, and Bionano Access Servers before server shipment and installation. IrysSolve servers, which have been in use with Bionano Access, also have Bionano Solve installed during an upgrade.

Bionano Solve (folder “tools”) is located at the /home/bionano directory on the compute server. The folder contains a collection of tools and scripts. Each individual tool is versioned independently. These tools together perform bioinformatics analyses on the compute server.

Compatibility

Bionano Solve 3.7 is compatible with Bionano Access 1.7 only.

Use of IrysSolve servers

If you are using an IrysSolve server to perform computational analysis on data generated by the Saphyr instrument, contact Bionano Support to reconfigure the IrysSolve prior to running any samples labeled with the DLS chemistry.

References

Visit <https://bionanogenomics.com/support-page/data-analysis-documentation/> for file format specifications and Theory of Operation documents.

Bionano EnFocus™ Fragile X Analysis Pipeline

- Initial release of targeted analysis of FMR1 repeat expansions that are related to Fragile X syndrome. Detailed performance specification and methods overview are described in the Bionano Solve Theory of Operation: Bionano EnFocus Fragile X Analysis Pipeline (PN 30457)

Bionano EnFocus™ FSHD Analysis Pipeline

- Fixed issue with identifying maps relevant to D4Z4 region that caused negative repeat counts to be reported
- Added ICS version to FSHD report
- Removed non-informative sections from report
- Added sample-level QC reporting based on molecule quality metrics and assessment of stable regions
- Updated report to clarify detection of CNVs overlapping SMCHD1
- Updated analysis to remove non-deterministic steps

AOH/LOH Detection

- Added detection of absence of heterozygosity (AOH)/loss of heterozygosity (LOH) events based on called zygosity of structural variants by the *de novo* assembly pipeline (not available as part of the Rare Variant Pipeline)
- Tested performance for detection of AOH/LOH > 15 Mbp in constitutional samples. Acceptable performance for > 25 Mbp was found and default filter was set accordingly.
- Implemented per variant estimate of probability that variant belongs to an AOH/LOH region
- Implemented initial model for confidence scoring of AOH/LOH regions based on size
- Added AOH/LOH stats to informatics report

Variant Allele Fraction estimation

- Implemented estimate of variant allele fraction for structural variants detected by *de novo* assembly and Rare Variant pipelines
- Implemented circular binary segmentation algorithm to detect segments of VAF similarity to aid visualization of aneuploidy and trisomy.

Copy number analysis pipeline

- Improved sensitivity of whole chromosomal aneuploidy detection

- Improved correction for systematic biases during coverage normalization
- Improved low VAF variant detection

De novo assembly pipeline

- Initial support for calling terminal deletions
- Preserve CNV and SV mask .bed files used in analysis in output .zip file
- Updated informatics reports to provide JSON version
- Add options to disable trimming of unlabeled molecule ends

Hybrid scaffolding

- Added option to override default N-base gap size

Molecule Quality Report (MQR)

- Standardized RefAligner parameterization to address differences in effective coverage estimates between MQR and assembly report
- Implemented MQR as JSON report to aid integration with external systems

Rare Variant Pipeline (RVP)

- Preserve CNV and SV mask .bed files used in analysis in output .zip file
- Updated informatics reports to provide JSON version to aid integration with external systems
- Added option to keep _full.xmap output for troubleshooting purposes

SV confidence

- Improved models for scoring translocation and inversion breakpoint confidence for human and non-human datasets using hyperparameter tuning and orthogonal SV data.
- Updated recommended confidence score filters based on new model

Variant Annotation Pipeline (VAP)

- Added annotation of CNV variant types
- Incorporated cross-reference of CNV events in duo and trio analyses.
- Added International System for Human Cytogenetic Nomenclature (ISCN) notation for all variant types (not aneuploidy, AOH, or triploidy)
- Updated cross-reference with control sample database to report frequencies by zygosity and ethnicity
- Recomputed DLE-1 control datasets for *de novo* assembly pipeline and RVP using Solve 3.7
- Streamline output by removing intermediate files from final .zip file

VCF conversion

- Added VAP annotation information to VCF for CNVs and SVs
- Standardized encoding of structural variant types and symbolic ALT alleles to meet VCF 4.2 spec
- Added LowConfidence and Mask FILTER values to synchronize with recommended confidence filtering levels in Access.
- Added INFO fields to preserve original Bionano structural variant type
- Added fractional copy number value for CNV gain and loss events
- Removed duplicate entries from SV calls on multiple maps
- Update SV size and orientation to synchronize with values reported in SMAP
- Added sample sex metadata fields
- Implement hemizygous genotype description on sex chromosomes
- Improve specification of breakpoint uncertainty
- Updated ALT allele of breakends to define location of mate per VCF spec

Miscellaneous

- Updated python2 software to python3
- Intra-chromosomal translocations renamed intra-chromosomal fusions
- Updated informatics reports to use JSON for improved chain-of-custody support and integration with external LIMS

Other known issues and limitations

- VAF calculation is performed independently from SV zygosity and AOH/LOH detection. We have observed cases where the LOH algorithm results and the VAF segments do not agree. This may occur because the analyses are performed independently using different rules for which variants to include in the calculation.

Solve 3.7.1 Updates

The following issues were addressed in the Solve 3.7.1 update:

Ticket	Description
BIOIN-2252	Addressed increase in inversion and duplication calls by improving SV calls formerly categorized as 'complex'
BIOIN-2343	Fixed bug in handling unexpected match group configuration in complex inversions during confidence scoring
BIOIN-2359	Removed hybrid scaffold dependency on specific Perl versions
BIOIN-2390	Fixed bug where insertions at the same position but with different lengths were filtered from VCF
BIOIN-2391	Fixed bug in confidence scoring and sizing of inversions
BIOIN-2420	Fixed missing dependency in runCharacterize.py
BIOIN-2424	Fixed bug in MQR by updating memory allocation to manage .bnx files with more than two billion molecules
BIOIN-2425	Added variant allele frequency to INFO field in VCF export
BIOIN-2428	Mask file dependency needed for custom mask creation
BIOIN-2437	Fixed bug in RefAligner assertion of obsolete condition which caused failures in alignmolvref stage for 1.5 Tbp human RVP
BIOIN-2474	Fixed bug where translocations were incorrectly flagged as _repeat
BIOIN-2475	Fixed bug in de novo pipeline in handling cleanup command line options (-c) so that full .xmap file not kept when requested
BIOIN-2480	Fixed RefAligner assertion error in handling label interval sample size of zero
BIOIN-2531	Fixed bug in de novo assembly without reference where configuration warnings caused pipeline failure
BIOIN-2556	Fragile X report missing Solve version

Technical Assistance

For technical assistance, contact Bionano Genomics Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

Type	Contact
Email	support@bionanogenomics.com
Phone	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7663
Website	www.bionanogenomics.com/support

Bionano Genomics, Inc.
9540 Towne Centre Drive, Suite 100
San Diego, CA 92121