



Bionano Genomics, Inc.
9540 Towne Centre Drive, Suite 100
San Diego, CA 92121
Phone +1 858 888 7600
Fax +1 858 888 7601
Support@bionanogenomics.com

Product Name: Bionano Prep Sp-G2 Blood & Cell Culture DNA Isolation Kit

CoA Number: CoA-00032

Part Number: 80060

Lot Number: 230405-0001885

Date of Manufacture (DOM): 3/21/2023

Expiration Date: 2/10/2024

Kit Components

Component Part Number	Component Name	Component Lot Numbers	Quantity
20442	RBC Lysis Buffer, 18 mL	221122-0000853	1 ea
20374	Cell Buffer, 50 mL	230316-0001886	1 ea
20443	Digestion Enhancer, 4 mL	230317-0001889	1 ea
20444	Lysis and Binding Buffer (LBB), 1.2 mL	230320-0001891	1 ea
20445	Wash Buffer 1 (WB1), 4.5 mL	230320-0001892	1 ea
20446	Wash Buffer 2 (WB2), 6.0 mL	230324-0001902	1 ea
20378	Elution Buffer (EB), 1.1 mL	230323-0001899	1 ea
20447	DE Detergent, 55 µL	230321-0001893	1 ea
20448	4mm Nanobind Disks, 12 ea.	230321-0001894	1 ea
20449	Protein LoBind Microcentrifuge Tubes, 1.5 mL, 12 ea.	230322-0001896	2 ea
20450	Protein LoBind Microcentrifuge Tubes, 0.5 mL, 12ea.	230322-0001897	1 ea
20451	SP Sheaths, 12 ea.	230323-0001900	1 ea
20452	Microcentrifuge Tubes, 2 mL, 12 ea.	230322-0001898	1 ea
20423	DNA Stabilizer, 350µL	230317-0001888	1 ea
20455	RNase A, 150 µL	230323-0001901	1 ea
20355	Ultrapure Water, 900 µL	230317-0001890	2 ea
20441	Thermolabile Proteinase K, 150 µL	230321-0001895	1 ea


QC Testing

Tests	Characteristic	Results
Functional Test	DNA recovered ≥ 39 ng/µL	☒ Pass ☐ Fail
	N50 (> 20 kbp) ≥ 150 kbp	☒ Pass ☐ Fail
	Map Rate $\geq 70\%$	☒ Pass ☐ Fail
	Label Density = 14 – 17 labels/100 kbp	☒ Pass ☐ Fail
	NLV $\leq 15\%$	☒ Pass ☐ Fail
	PLV $\leq 10\%$	☒ Pass ☐ Fail



Bionano Genomics, Inc.
9540 Towne Centre Drive, Suite 100
San Diego, CA 92121
Phone +1 858 888 7600
Fax +1 858 888 7601
Support@bionanogenomics.com

Bionano Genomics hereby certifies that the quality of this product conforms with internal specifications and meets required quality standards. This is to certify that the information listed above is a true and accurate copy of the data obtained for this material. This product is intended for research use only and not for use in diagnostic procedures.

Bionano Genomics Quality Assurance	
Review and Approved By:	
	3-23-23
Signature	Date