



# **Bionano Image Detect™ v1.0**

## **Release Notes**

Document Number: 30162

Document Revision: A

# Table of Contents

Introduction .....	3
Known Issues and Limitations .....	3
Technical Assistance.....	4

## Introduction

This document describes the release of Bionano Image Detect™ v1.0.

This is an initial release.

## Known Issues and Limitations

- Bionano Image Detect™ is only compatible with the Saphyr system.
- Only green-labeled data is supported.

## Technical Assistance

For technical assistance, contact Bionano Genomics Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support page or by request through e-mail and telephone.

Type	Contact
Email	support@bionanogenomics.com
Phone	<b>Hours of Operation:</b> Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7600
Website	<a href="http://www.bionanogenomics.com/support">www.bionanogenomics.com/support</a>